

Phonak Audéo™ E.

Product information

Audéo E is the next generation Receiver-In-Canal (RIC) portfolio designed for mild to profound hearing losses. Audéo E is built for the moments that matter most, designed for adaptive clarity and seamless connections to the voices and moments that shape your world.



Product description



Audéo E-Sphere
Audéo E-Sphere Trial
E90/E70/E50

Audéo E-R
Audéo E-R Trial
E90/E70/E50/E30

Battery		built-in rechargeable
Multi-function button		•
Bluetooth®		•
Hydrophobic coating		•
IP rating		IP68 ¹
Dimensions	L x W x D	28.7 x 12.6 x 8.9 (1.13 x 0.5 x 0.35")
	Weight	2.4 g (0.085 oz)
Charging times ²	0% to 100%	2.5 h
	0% to 90%	2 h
	0% to 55%	1 h
	30 min. fitting	15 min.
S Receiver 6.0		
Max. Power Output (dB SPL)	2 cc coupler	112
Max. gain (dB)	2 cc coupler	47
Frequency range (Hz)		<100 - 10000
M Receiver 6.0		
Maximum Power Output (dB SPL)	2 cc coupler	114
Max. gain (dB)	2 cc coupler	51
Frequency range (Hz)		<100 - 10000
MAV Receiver 6.0		
Maximum Power Output (dB SPL)	2 cc coupler	119
Max. gain (dB)	2 cc coupler	50
Frequency range (Hz)		<100 - 9000
P Receiver 6.0		
Maximum Power Output (dB SPL)	2 cc coupler	122
Max. gain (dB)	2 cc coupler	59
Frequency range (Hz)		<100 - 7500
UP Receiver 6.0		
Max. Power Output (dB SPL)	2 cc coupler	132
Max. gain (dB)	2 cc coupler	71
Frequency range (Hz)		<100 - 8000

Housing colors



P1
Sand
Beige



P3
Sandalwood



P4
Chestnut



P5
Champagne



P6
Silver
Gray



P7
Graphite
Gray



P8
Velvet
Black



8A
Sienna
Red

¹ IP68 indicates that the hearing aid is water and dust resistant. It survived continuous immersion in 1 meter of fresh water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

² Actual charging time depends on the remaining battery life.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

Performance levels

	Audéo E-Sphere			Audéo E-R			
	Premium	Advanced	Standard	Premium	Advanced	Standard	Essential
AutoSense OS™ AI 8.0							
Spheric Speech in Loud Noise	•	•					
Speech in Loud Noise				•			
Music	•	•		•	•		
Speech in Car	•	•	•	•	•	•	
Comfort in Noise	•	•	•	•	•	•	
Speech in Noise	•	•	•	•	•	•	•
Calm Situation	•	•	•	•	•	•	•
Media speech + mic	•	•	•	•	•	•	•
Media music + mic	•	•	•	•	•	•	•
Additional programs							
Max. additional programs	4	4	4	4	4	4	4
Spheric Speech in Loud Noise	•	•	•				
Speech in Loud Noise	•	•		•	•		
Speech in 360°	•	•		•	•		
Speech in Car	•	•	•	•	•	•	
Speech in Noise	•	•	•	•	•	•	•
Calm Situation	•	•	•	•	•	•	•
Comfort in Noise	•	•	•	•	•	•	•
Music	•	•	•	•	•	•	•
Acoustic phone	•	•	•	•	•	•	•
Custom program	•	•	•	•	•	•	•
Streaming programs							
RogerDirect + mic	•	•	•	•	•	•	•
PartnerMic + mic	•	•	•	•	•	•	•
Phone call + mic	•	•	•	•	•	•	•
Features							
Spheric Speech Clarity 3.0	•	•	•				
Speech Enhancer	•			•			
SpeechSensor	•	•		•			
StereoZoom 2.0	•	•		•	•		
Dynamic Noise Cancellation	•	•		•	•		
Motion Sensor Hearing	•	•	•	•	•	•	
RogerDirect™	•	•	•	•	•	•	•
ActiveVent™	•	•	•	•	•	•	•
Disable Hands-free Calls¹	•	•	•	•	•	•	•
Tap control							
- Siri/Google Assistant	•	•	•	•	•	•	•
- Pause/resume streaming	•	•	•	•	•	•	•
- Accept/end phone call	•	•	•	•	•	•	•
SoundRelax	•	•	•	•	•	•	
DuoPhone	•	•	•	•	•	•	•
Real Ear Sound	•	•	•	•	•	•	•
UltraZoom	•	•	•	•	•	•	•
Health functionalities²	•	•	•	•	•	•	•
WindBlock 2.0	•	•	•	•	•	•	•
User Preference Tuning	•	•	•	•	•	•	•
Tinnitus Balance	•	•	•	•	•	•	•
Auto Acclimatization	•	•	•	•	•	•	•
Environmental Balance	•	•	•	•	•	•	•
SoundRecover2	•	•	•	•	•	•	•
WhistleBlock	•	•	•	•	•	•	•
NoiseBlock	•	•	•	•	•	•	•
QuickSync	•	•	•	•	•	•	•
AOV	•	•	•	•	•	•	•
Fine tuning channels	20	20	16	20	20	16	12

¹ For high security environments, such as SCIF
² Available in combination with myPhonak

CROS compatibility

	Audéo E-Sphere Audéo E-Sphere Trial	Audéo E-R Audéo E-R Trial
Phonak CROS E-R		•


Fitting

	Audéo E-Sphere Audéo E-Sphere Trial	Audéo E-R Audéo E-R Trial
Phonak Target	12.0 or higher	12.0 or higher
Noahlink Wireless, Noahlink Wireless 2	•	•

Acoustic coupling

Length	S Receiver	M Receiver	MAV Receiver	P Receiver	UP Receiver
Available in five different lengths (00, 0, 1, 2, 3), L and R	•	•	• ¹	•	• ¹
Dome					
EasyGuard™ Open*	•	•		•	
Available in three sizes (S, M, L)					
EasyGuard™ Vented*	•	•		•	
Available in three sizes (S, M, L)					
EasyGuard™ Power*	•	•		•	
Available in five sizes (S, S/M, M, M/L, L)					
Cap dome	•	•		•	
SlimTip					
Acrylic	•	•		•	
Silicone	•	•		•	
Titanium AV			• ¹		
Universal AV Available in four shapes (A, B, C, D)			• ¹		
cShell					
Acrylic	• ¹	• ¹		• ¹	• ¹
Titanium	• ¹	• ¹		• ¹	• ¹

Fitting ranges



Mild to moderately-severe hearing loss.

Mild to severe hearing loss.

Mild to severe hearing loss.

Mild to severe hearing loss.

Moderate to profound hearing loss.

¹ Length option 00 not available

* Also available in standard (previous generation) domes in three sizes (S, M, L)

Phonak Audéo E charging options



Phonak Charger RIC E



Phonak ChargerGo RIC I E

Compatible with Audéo E-Sphere, Audéo E-Sphere Trial Audéo E-R, Audéo E-R Trial CROS E-R	•	•
Power indication LED	•	•
Battery indication LED		•

Phonak Audéo E wireless technology

Direct calls/audio streaming	Bluetooth® Classic Audio
Auracast™ broadcast	Bluetooth® LE Audio

Wireless communication portfolio

Audéo E-Sphere
Audéo E-Sphere Trial
Audéo E-R
Audéo E-R Trial


Wireless accessory	
TV Connector	•
Phonak PartnerMic	•
Phonak RemoteControl	•
Roger	
Roger microphone portfolio	•
Apps	
myPhonak	•

For more information, visit www.phonak.com/audéo-e

For information on Roger microphones, visit www.phonak.com/roger

For more information on myPhonak and myPhonak Junior, visit www.phonak.com/apps

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by Bluetooth SIG, Inc. Any use of such marks by Sonova AG is under license.



life is on

Since 1947, Phonak is dedicated to preserving life quality by opening new acoustic worlds. Already back then, in the days of our foundation, our company was driven by a formative conviction: We believe that well-hearing equates to well-being and thus is essential for living life to the fullest. In fact, the sense of hearing is directly linked to social, emotional, cognitive and physical well-being. Today as in future, we strive to offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

phonak.com



027-0742-02/V2.00/2026-07/jb © 2026 Sonova AG. All rights reserved.