

For professionals

ReSound GN

ReSound Sensia™

This is what the
best sounds like



GN Bringing People Closer

A revolution in hearing technology

Powered by our new AI platform, ReSound Sensia introduces Intelligent All-Around*, an automatic feature that uses a powerful trio of technologies to deliver the best hearing clarity in both noise and quiet.



AcoustIQ™
Most accurate real-time sound sensing technology¹



Clear focus
Narrowest 4-mic beamforming⁷



DNN Denoising
17x more efficient²

Unlock the full sound environment

AcoustIQ senses and adapts to the sound environment, engaging the right features at the right time.

Superior access to all-around sound, preferred 90% of the time by users.³

88% greater accuracy in real-time classification and feature engagement than the closest competitor.¹

91% accuracy in classifying the environment.⁴

Focus on speech – instantly

Clear Focus creates a forward-facing beam, so users can hear what they want simply by looking toward it.

Preferred by users 93% of the time.⁵

13% SNR improvement compared to legacy.⁶

83% of users experience superior speech clarity in noise compared to the closest competitor.⁷

Follow conversations with less effort

A new approach to **DNN denoising** that steps in automatically to keep conversations clear and make background noise feel less intrusive.

Validated by Hörzentrum Oldenburg:

ReSound Sensia is the only solution that combines 4-mic beamforming and DNN noise reduction simultaneously for better audibility and lower listening effort in noise.⁸⁻⁹

Resulting in best audibility of the intended speaker for 80% of users.¹⁰

From our portfolio to yours

Bring what the best sounds like to all your clients with a new family of Premium hearing aids built on years of innovation so your clients can hear exactly what they want to hear.



Next-level clarity on all levels

ReSound Sensia brings an externally validated, industry-leading sound experience to more users by extending our latest AI-driven technology across programs and tech levels, so more of your clients can benefit from clearer conversations in challenging listening environments.



The smallest AI rechargeable RIE

From the world's smallest rechargeable AI RIE¹¹ to a full range of Customs, a CROS solution, and Bimodal connectivity, you can offer your clients a line-up of IP68-rated hearing aids that bridge all-day comfort with all-weather durability.



Ready, set, connect with Auracast™ broadcast audio

Bring the clearest streaming experience¹² with Bluetooth® Low Energy (LE) Audio and Auracast™ to your clients, connecting them to a solution preferred by 87% of users compared with telecoil.¹³



Find confidence in ReSound Sensia

ReSound hearing aids are recognised not only by leading institutions but by the people who rely on them every day.

Tested.

Rigorously tested in real-world scenarios.

Preferred.

Tried and preferred over competition by users.

Proven.

Validated by leading, world-renowned institutions for hearing in noise.

Learn more



¹Groth & Cui (2026). ²GN Proprietary data on file (2026). ³Jespersen & Skovlund (2025). ⁴Groth & Cui (2024). ⁵Jespersen & Skovlund (2025). ⁶Jespersen & Hashem (2026). ⁷Jespersen, Dieu & Rubachandran (2024). ⁸Jespersen, Heeren & Quilter (2026). ⁹Quilter, Heeren & Jespersen (2026). ¹⁰Compared to other premium hearing aid brand with DNN denoising. Jespersen & Rubachandran (2026). ¹¹GN Proprietary data on file (2026). ¹²Bluetooth SIG. ¹³Jespersen, Hashem, Aranda de Toro & Fernandez (2026).

ReSound Sensia style options


ReSound GN



microRIE 60S
Rechargeable

RIE 61
312 battery


RIE 62
13 battery



CIC-NW
10A battery

ITC
312 battery

ITC/ITE
Rechargeable



Key feature overview

For a detailed feature overview across styles, please refer to the data sheets

Features	9	7	5	4
Speech Intelligibility & Awareness	★★★★★	★★★★	★★★	★★
All-Around Directionality	Intelligent All-Around* AI+	360 All-Around	All Access Directionality	Binaural Directionality III
Speech Intelligibility in Noise	Intelligent Focus*/Clear Focus AI+		Front Focus	Ultra Focus
Comfort	★★★★★	★★★★	★★★	★★
DNN Noise Reduction (Intelligent Noise Tracker)	Mild – Strong* AI+	On/Off* AI+		
Noise Reduction (Noise Tracker II)	Per environment (-10dB)	Strong (-10dB)	Considerable (-8dB)	Moderate (-6dB)
Automatic Volume Adjustments	Environmental Optimizer II	Environmental Optimizer		
Wind Noise Reduction	Mild - Strong		Mild - Moderate	
Sound Processing	★★★★★	★★★★	★★★	★★
Environmental Classification & Feature Engagement System	AcoustiQ AI+			
Warp Bands	17	14	12	
Feedback Management	DFS Ultra III AI+			
CROS/BICROS	Yes*			
Bimodal Compatibility	Yes		Yes	
Connectivity ★★★★★				
 MFi & Direct Audio Streaming with Android  Hands-free Calls**  Bluetooth LE Audio & Auracast***  ReSound Smart 3D™ App  Streaming Bass Boost  Multi-Mic+ & TV-Streamer+				
In-clinic ecosystem ★★★★★				
 ReSound Smart Fit™ 2.4  Flexible Demo Devices  Personalisable Patient Experience Levels  NAL-NL3™ & Audiogram+  4 Customizable Programs  AllFit Domes for better acoustic fit				

* Only available for the rechargeable microRIE style.

** Not available for CIC and battery-type ITC styles. Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation), iPad (10th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later. Bluetooth® Low Energy (LE) Audio works with ReSound Sensia and compatible mobile devices with Bluetooth 5.3 or later. See compatibility on gnhearing.com/products/compatibility.

*** Not available for CIC style.

Worldwide Headquarters

GN Hearing A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: +45 4575 1111
gnhearing.com
CVR no. 55082715

United Kingdom

GN Hearing UK Ltd.
Unit 13
Talisman Business Centre
Bicester
Oxon OX26 6HR
Tel: +44 1869 352800
gnhearing.com

Australia

GN Hearing Australia Pty Ltd
Suite 1001, Level 1
1 Julius Avenue
North Ryde, NSW 2113
Australia
Tel.: (free) 1800 658 955
gnhearing.com

New Zealand

GN Hearing New Zealand
Limited
Ground Floor, North Entrance
4 Fred Thomas Drive
Takapuna, Auckland, 0622
Tel.: (free) 0800 900 126
gnhearing.com



© 2026 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone and iPad are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast word mark and logos are trademarks owned by Bluetooth SIG, Inc. Cochlear, Hear now. And always, Nucleus, Kanso, Baha, the elliptical logo, and marks bearing an ® or ™ symbol, are either trademarks or registered trademarks of Cochlear Limited or Cochlear Bone Anchored Solutions AB (unless otherwise noted).